

OWNER'S MANUAL



THE "GO-ANYWHERE" SPORT

"SPT"



US Patent No. 29/389,689 (ref. website for international patents)

IMPORTANT - READ THIS CAREFULLY AND COMPLETELY!

Congratulations and Welcome to GO! Mobility Solutions' family of innovative products and aids designed to assist physically-challenged individuals to get out and GO! in ways never thought of before. We are confident that this product will change your life as it has for so many others. Please be advised that, from time-to-time, improvements and modifications are incorporated into our products which may not be reflected in the printed manual. For the most up-to-date information please logon to www.GoesAnywhere.com, click on the "Products" tab and find the downloadable link to the User Manual PDF file for the product that you own. Should you have any questions or require any additional information not found within this manual, please logon to the website and click on the "Contact Us" link or phone us, toll free, at (800) 359-4021. ***It is highly recommended that you review this manual carefully prior to use of your GO-Anywhere Chair™.*** We would also appreciate receiving any feedback and comments that you may have to offer so that we may continue to improve our products and our customers' quality of life.

WARRANTY NOTICE – WARNING AND DISCLAIMER

Please see LIMITED WARRANTY below. If this warranty is not acceptable, return the chair for a full refund PRIOR TO USE. Contact GO! Mobility Solutions for shipping instructions and a Return Authorization BEFORE returning any product.

The "GO-Anywhere" Chair SPT™ has been designed and successfully tested for use by persons NOT exceeding 300 pounds. This product is NOT intended for use by persons weighing more than 300 pounds. Should any individual exceeding 300 pounds choose to use this product for any reason, he or she does so entirely at his or her own risk and contradictory to the recommendations and warnings of GO! Mobility Solutions. Use

by any person in excess of the 300-pound weight limit shall void this warranty in its entirety.

The “GO-Anywhere” Chair™ may become damaged in its soft travel case should any heavy object be dropped on it.

If any aluminum tubing becomes dented, telescoping parts may be affected rendering them dysfunctional. These parts should be replaced as soon as possible.

Carefully follow the step-by-step unpacking, assembly, operating and disassembly instructions as specified below. For additional assistance please view the corresponding videos at www.GoesAnywhere.com.

PRODUCT OVERVIEW


The GO-Anywhere Sport™ is our newest addition to GO! Mobility Solutions’ portable GO-Anywhere Chair™ technology offering a new bathing solution for active, independent wheelchair users. The Sport has been designed on a versatile frame so that it may be used in a three different ways:

In its most compact configuration it serves as a lightweight, but strong Tub Bench, open to both sides for easy access beneath, comes with height-adjustable legs and a stylish and convenient carrying case.

A separate (optional) Transfer Section accessory utilizes two of the legs from the Tub Bench, which are relocated from the inside of the tub to the outside. It has been thoughtfully designed to work in both left- and right-handed tubs. The Transfer Section accessory comes with its own travel bag which attaches very simply to the main travel bag.

An (optional) Backrest is offered as another accessory for those requiring lumbar support. The backrest fits snugly into the carrying case provided along with the tub bench.

Assembly and disassembly are a snap, require no tools and can be completed in a matter of seconds. The seat height may be adjusted by shortening or extending the telescoping legs using spring pin buttons and a series of adjustment holes.



Weight-bearing capacity = 300 lbs maximum.

Weight (in bag) = 19 lbs + 9 lbs (Transfer Accessory in bag)

Bag Dims: Tub Chair: 7.5"x 20"x 11" / Transfer: 5.5"x 20"x 11"

The GO-Anywhere Sport™ is:

- Safe and easy to use.
- Simple to assemble; no tools are required.
- Lightweight and durable, constructed of high-quality anodized aircraft-grade aluminum.

- Compact, Convenient and PORTABLE! Meets U.S. FAA criteria as “carry-on” baggage (<45 lineal inches).
- Rust and corrosion-resistant, easy to clean and virtually maintenance-free.
- Adjustable; the seat height may be easily adjusted without the use of any tools.
- Comfortable; the seat and backrest cushions are constructed of durable, waterproof foam, providing for maximum comfort. Both are attached to the chair frame with Velcro™ rendering them easily removable for cleaning or replacement.

PRACTICE MAKES PERFECT!

PLEASE do not become discouraged if you have problems either with the assembly or use of the “GO-Anywhere” Chair™ on your first try. Like with anything else, familiarity and practice will serve to continually improve your experience with your new chair. Additionally, as your chair is brand new, some of the hinges or telescoping parts will be snug, initially. They will “break-in” and operate more easily with each use. Of course, should you have ANY problems, please log on to www.goesanywhere.com, e-mail us at support@goesanywhere.com or call us at (800) 359-4021.

UNPACKING INSTRUCTIONS AND PARTS INVENTORY



1. *PAY CLOSE ATTENTION TO THE ORDER IN WHICH YOU UNPACK THE CONTENTS OF THE CASE SO THAT YOU WILL BE ABLE TO REPACK IT EASILY IN THE REVERSE ORDER.*

2. The “GO-Anywhere” Sport™ is composed of twelve (12) individual components in addition to the (optional) Transfer Section and Backrest accessories:

(1) **travel bag**

(1) folding **tub bench**

(2) **seat cushions** (pre-assembled w/ Tub Bench)

(4) (telescoping) **upper legs**

(4) (telescoping) **lower legs w/ adj. non-skid glides**

Please note that, although they are delivered as one pre-assembled unit, the **tub bench** and **seat cushions** are separate parts. Similarly, the four **upper legs** and four **lower legs w/ adj. non-skid glides** may also delivered pre-assembled. *EVEN THOUGH YOU MAY HAVE AS FEW AS FIVE (5) PARTS TO ASSEMBLE, PLEASE MAKE SURE THAT ALL TWELVE (12) COMPONENTS ARE ACCOUNTED FOR.*

ASSEMBLY AND USE

1. Unzip and open the lid of the travel bag and remove the folded **tub bench** and four **telescoping legs**. (If you purchased the (optional) **transfer section** you will have received two travel bags. The **tub bench** is in the larger, wheeled bag; the **transfer section** is in the smaller shoulder bag.)

2. Unfold the seat and install the through pin which is attached to the frame via nylon lanyard.

3. Insert each of the four **lower legs w/ adj. non-skid glides** into the **upper legs** using the spring pin buttons in the lowers.

4. If using the **tub bench WITHOUT** the (optional) **transport section**, install each of the four assembled **telescoping legs** into their respective leg sockets on the underside of the frame, giving each leg a firm push until you feel it snap into place. The tension of the ball plunger(s) located at each of the leg sockets may be loosened or tightened as needed; OR

5. If attaching the **tub bench to the (optional) transfer section**, install only two **telescoping legs** into the **tub frame** (using the leg sockets in the frame that will be placed closet to the wall).

Place the **tub bench** into the bathtub with the installed legs at the far end of the tub and allow the side of the **tub bench** having the unused leg sockets to rest on top of the outside edge of the tub.



After removing the **transfer section** from its bag, install the two remaining **telescoping legs** into their respective sockets on the frame. Turn the

transfer section over and insert the “hooks” beneath the **tub section** frame. Push the **transfer section** into the **tub section**.

6. Once they are joined, install the through pin which is attached to the frame via nylon lanyard, ensuring that it



penetrates both tub bench and transfer section frames.

7. If the (optional) **backrest** was purchased, it may now be installed.



Using the adjustable non-skid swivel glides at the bottom of each leg, adjust the glides so that they raise or lower until the chair is stable within the tub and does not wobble or rock. Loosen the hex nut on each of the glides first, then adjust the length of the bolts and finally re-tighten the hex nut to lock the length. Ensure that the bench is secure before transferring to it.



Note: The swivel feet have a 1+” adjustment range. Make big adjustments with the spring pins and fine adjustments with the swivel feet. ALWAYS CHECK THE SWIVEL FEET BEFORE USING THE CHAIR. Extending them too far out of their leg might cause them to be unstable or, even worse, to dislodge from the leg.

DISASSEMBLY & PACKING INSTRUCTIONS

1. Remove the through pin which secures the **tub bench** and **transfer section** frames. Disconnect the frames.
2. Remove all four legs from their respective sockets by giving each a firm tug. Depress the spring pin buttons and slide the smaller-diameter **lower legs** as far in as they will go, into the wider-diameter **upper legs**, making them as short as possible.

3. Fold the **tub bench** frame into its closed position and slide it into the travel bag with its hinged-side at the bottom. Place all four shortened leg assemblies in between the folded halves and zip the **travel bag** closed.

NOTE: When packing the (optional) **backrest**, REMOVE the backrest cushion and place the backrest (with the frame supports down) into the **travel bag** after packing the folded **tub bench**. Separate the four **upper legs** from the four **lower legs**, place the **upper legs** into the bag vertically, followed by the **lower legs**, horizontally. Place the **backrest** cushion on top and zip the travel bag closed.



MAINTENANCE INSTRUCTIONS

1. Following each use, towel-dry all parts of the chair. If used as an everyday-use chair, it is recommended that after each use, your GO-Anywhere Chair™ be laid on its side for a little while to permit any accumulated moisture to drain out.

2. From time-to-time, carefully remove the seat cushions, so as not to tear the front face of the cushions from the seat. Dry both the aluminum base and the backs of the cushions to prevent the growth of any mold or mildew. Make sure to peel the cushions off only when they are dry to prevent having the Velcro peel off with them.

Hydrogen peroxide sprayed on the cushions works very well to clean them or you may wipe them down with an antibacterial soap.

3. Replacement cushions may be ordered from GO! Mobility Solutions either on-line at www.goesanywhere.com or by calling (800) 359-4021.

NOTIFICATION OF PATENTS

Please refer to www.GoesAnywhere.com to review all US and international patents and trademarks.

TRAVEL WITH YOUR GO-ANYWHERE CHAIR™

When traveling entails flying and airport security, be prepared! You are, most assuredly, going to be confronted with one or more potential obstacles.

Know your rights! There are numerous resources available on-line from which you may easily glean the specifics using a keyword search for “airline and wheelchair”. Your GO-Anywhere Chair™ is a “fragile medical device” which should be “gate-checked” *whenever possible* and is not chargeable as passenger baggage.

You can review a chart of all airlines baggage fees at:

http://www.bestfares.com/news-newsID589-Airline_Checked_Baggage_Fee_Chart.html.

Review TSA regulations specific to travelers with disabilities at:

<http://www.tsa.gov/travelers/airtravel/specialneeds/index.shtm>.

Air travel for people with disabilities is a hassle and can be embarrassing and downright uncomfortable. Following the steps in the aforementioned process will not change this fact, unfortunately, though it will make it a little less miserable.

And, because your rights are protected by the Air Carrier Access Act (14 CFR, Pt 382) if you experience any sort of problem that is directed solely to your disability, immediately ask for a Complaint Resolution Official. The C.R.O. is there to resolve problems that may arise due to your disability and for no other purpose. They are, by law, supposed to be available to any departing or arriving flight. Do not be forced into any "assistance" that you did not ask for or do not want and on the other hand, ask for assistance that is supposed to be there but is not readily available. Go to the following link for more information:

http://www.disabilitytravel.com/airlines/air_carrier_act_details.htm

FOR ADDITIONAL UP-TO-DATE INFORMATION AND INSTRUCTIONS PLEASE VISIT OUR WEBSITE AT:

www.GoesAnywhere.com.

LIMITED WARRANTY AND RETURN POLICIES

Please read this BEFORE operating or using your GO-Anywhere Chair™.

By operating or using the GO-Anywhere Chair, you agree to the terms of this warranty. If you do not agree, return the chair to GO! Mobility Solutions BEFORE use to receive a full refund of the price you paid. Thereafter, this product may be returned at any time within 30 days following delivery for a full refund, less a 20% restocking fee and at your own shipping expense (assuming the product is in "like new" condition), and less a \$50 non-refundable fee for the cushions which are yours to keep. Health codes prohibit the resale of cushions, used or not.

Should a "Trial-Use Period" be advertised and offered it shall modify the previous paragraph accordingly.

How long is the limited warranty period?

GO! Mobility Solutions warrants its products to be free from defects in material and workmanship for up to five years, frame; twelve months, non-frame parts and one year labor starting from the date of original purchase from GO! Mobility. This warranty excludes foam cushions beyond thirty days. This limited warranty also applies to repaired or replacement Products for a period of either (a) 90 days from the date the repaired or replacement Product is shipped to you or (b) the time remaining on the original limited warranty, whichever is longer.

What does GO! Mobility Solutions' limited warranty cover?

The manufacturer of this GO! Mobility Solutions product, GO! Mobility Solutions warrants to the holder of a valid proof of purchase ("Consumer" or "you") that the product and all accessories provided by GO! Mobility Solutions in sales package ("Product") are free from material defects in material and workmanship, pursuant to the following terms and conditions, when used normally and in accordance with operation instructions.

This limited warranty does not cover:

1. Product that has been subjected to misuse, accidental damage, shipping or other physical damage, abnormal operation or handling, neglect, fire, excessive heat, natural disaster; or
2. Product that has been damaged due to repair, alteration or modification by anyone other than an authorized service representative of GO! Mobility Solutions; or
3. Product whose warranty/quality stickers, Product serial numbers plates or electronic serial numbers have been removed, altered or rendered illegible; or
4. Product used for commercial or institutional purposes (including but not limited to Products used for rental purposes); or
5. Product returned without valid proof of purchase (see 2 below); or
6. Foam and vinyl (cushions), beyond a thirty (30) day period, due to their fragile nature; or
7. SOFT TRAVEL CASE – DAMAGE MAY OCCUR TO THE ALUMINUM CONTENTS IF DROPPED OR IF OBJECTS ARE DROPPED ON IT. DO NOT CHECK LUGGAGE AT AIRPORTS. GATE CHECK UNLESS EXPRESSLY FORBIDDEN.

What will GO! Mobility Solutions do if the Product is not free from material defects in materials and workmanship during the limited warranty period ("Materially Defective Product")?

If the product is received within the warranty period and if the product is found by to be defective within the terms of this warranty GO! Mobility Solutions will repair or replace the defective product at its option and cost.

How do you get warranty service?

1. To obtain warranty service in the United States of America. Call 1-800-359-4021 for instructions regarding where, and what, to return. Before calling for service, please check the user's manual.
2. Please do not return the chair to GO! Mobility Solutions without prior authorization from GO! Mobility Solutions. Any shipping and handling expense is the customer's responsibility unless otherwise authorized by GO! Mobility Solutions.
3. Except as provided by applicable law, you assume the risk of loss or damage during transit and transportation and are responsible for delivery or handling charges incurred in the transport of Product(s) to the service location. GO! Mobility Solutions will return repaired or replaced product under this limited warranty to you, transportation, delivery or handling charges prepaid. GO! Mobility Solutions assumes no risk for damage or loss of the Product in transit.
4. If the Product failure is not covered by this limited warranty, or proof of purchase does not meet the terms of this limited warranty, GO! Mobility Solutions will notify you and will request that you authorize the cost of repair and return shipping costs for the repair of Products that not covered by this limited warranty.

What must you return with the Product to get warranty service?

1. Return the product, in whole or in part, as was determined during the telephone call you made to report the issue requiring service to GO! Mobility Solutions.
2. Include a written description of the malfunction or difficulty;
3. Include "valid proof of purchase" (sales receipt) identifying the Product purchased (Product model) and the date of purchase or receipt; and
4. Provide your name, complete and correct mailing address, and telephone number.

Limitations and Exclusions: Implied warranties, including those of fitness for a particular purpose and merchantability (an unwritten warranty that the product is fit for ordinary use) are limited to one year from date of purchase. Some states do not allow limitations on how long an implied warranty lasts, so the above limitation may not apply. GO! Mobility Solutions is obligated to provide the purchaser with repair or replacement as described above so long as the warranty period is in force. This exclusive warranty remedy will not have failed as long as GO! Mobility Solutions is willing and able to repair or replace the product as described. Should this remedy be held to have failed, the only remaining warranty obligation of GO! Mobility Solutions shall be to provide the customer with a full refund. In no event shall GO! Mobility Solutions be liable for any indirect, special, incidental, consequential, or similar damages (including, but not limited to lost profits or revenue, inability to use the product, or other associated equipment, the cost of substitute equipment, and claims by third parties) resulting from the use of this product. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

This warranty is the complete and exclusive agreement between you and GO! Mobility Solutions. It supersedes all other written or oral communications related to this product. GO! Mobility Solutions provides no other warranties for this product. The warranty exclusively describes all of GO! Mobility Solutions' responsibilities regarding the product.

There are no other expressed warranties. No one is authorized to make modifications to this warranty and you should not rely on any such modification.

State Law Rights: This warranty gives you specific legal rights, and you also have other rights which vary from state to state.

DISCLAIMER

Except for the warranty as provided herein, and the acknowledgement by GO! Mobility Solutions that the chair, as manufactured by it, is fit for the general purpose for which

most persons acquire a chair of its kind, GO! Mobility Solutions provides that you accept the chair as is, without warranties, either express or implied. GO! Mobility Solutions makes no warranty of fitness for your particular purpose and no warranty of merchantability beyond that already stated. No warranties extend beyond the duration of the express warranty stated above.